

	Sonosite	SonixHub	Q-path	UltraLinq
wireless upload	+	+	+	+
local server	+	+	+	
cloud storage			+	+
transfer to PACS	+	+	+	
remote access	+	+	+	+
multi-vendor support	+	+	+	+
custom worksheets		+	+	
integrated worksheets		+ ^I	+ ^{II}	
populate patient info	+	+	+	+
send QA emails	+	+	+	+
export as PDF	+	+	+	+
export via HL7	+	+	+	+
billing reports	+	+	+	+
diagnostic vs educational		+	+	+
user statistics	+	+	+	+
Excel export option	+	+	+	+
teaching database		+	+	+
image de-identification	+	+	+	+ ^{III}
storage back-up	+	+	+	+

wireless upload	<i>wireless transfer of images from ultrasound machine to server</i>
local server	<i>scan images and reports are archived on a local server</i>
cloud storage	<i>scan images and reports are archived on a remote server</i>
transfer to PACS	<i>scan images can be transferred to the Radiology PACS system</i>
remote access	<i>scans can be viewed outside the hospital network (VPN access required)*</i>
multi-vendor support	<i>able to import scans from multiple machine vendors</i>
custom worksheets	<i>scan interpretation and QA worksheets can be customized</i>
integrated worksheets	<i>interpretations can be entered directly on ultrasound machines</i>
populate patient info	<i>patient information auto-populated into ultrasound machines</i>
send QA emails	<i>send QA emails with reports generated from completed worksheets</i>
export as PDF	<i>scan reports can be exported in PDF format</i>
export via HL7	<i>scan reports can be transferred to the EMR via HL7 standard format</i>
billing reports	<i>scan reports can be exported for billing purposes</i>
diagnostic vs educational	<i>separate purely educational from diagnostic scans within the database</i>
user statistics	<i>tracks sonographer scan numbers and accuracy statistics</i>
Excel export option	<i>scan statistics can be exported into standard Excel spreadsheets</i>
teaching database	<i>exemplary scans can be tagged within the database for teaching purposes</i>
image de-identification	<i>automated removal of patient identifiers when scan images are exported</i>
storage back-up	<i>perform periodic backups of the scan database</i>

Key

- * VPN required except Q-path Cloud and UltraLinq
- I available only on Ultrasonix ultrasound machines
- II available only on Zonare ultrasound machines
- III masks must be manually applied, not automated